

# Instrument Standard Operating Procedure

---

Raman Reflex

*Last Updated: 6/25/2020*

**(Under Construction)**

You must be a “Qualified Self-User” to operate this instrument independently.

You must be on the labs “Instrument Reservation Schedule” before touching the instrument for any reason.

Any problems, STOP, Post a note on the instrument and send an email to [mtim@mit.edu](mailto:mtim@mit.edu) immediately.

Do not perform any maintenance.

Do not adjust any optics.

## **Instrument Hazards:**

Lasers, Cryogenic Liquids Chemicals & Electricity.

## **Laser Safety:**

Anyone operating the laser must have attended the MIT Radiation Protection Offices Laser Safety Training. Call X2-EHSS to schedule a class.

Lasers: Class 3B, Interlocked system (No exposed laser beams).

Wavelength and Power: 785nm 100mw, 533nm 50mw, 437nm 25mw.

## **Required Apparel:**

Safety Glasses, Laser Safety Glasses when interlocks are disabled.

Wear clean gloves when handling optics.

## **Check out the system when you arrive:**

If you encounter any problems: Stop work and send Tim ([mtim@mit.edu](mailto:mtim@mit.edu)) an email immediately.

1. Lasers, Spectrometer and Microscope Laser Safety Enclosures – Closed.
2. Computer – On
3. Software - Closed.
4. Spectrometer Power – Always On

Operation:

**Logon to computer:**

User name: Admin

Password: 134139

Open the software - "Wire".

Turn the Lasers on (Key).

785nm laser must warmup ~60 seconds before performing health check.

**Verify Performance:**

Run the Health Check

On the machine status toolbar:

Choose your desired Laser Wavelength & Grating.

473nm & 532nm – Regular or Confocal

785nm – Line Focus Only

**Perform your work:**

**Save your data:**

Walk away with a copy of your data if it is important to you.

Use a USB or the computers internet connection to transfer your data.

**Data Analysis:**

Wire software is available to users.

Available on the desktop of the Raman computer.

## **Shutdown:**

Close the Software  
Turn off the Lasers (Key)

Cleanup.  
Disengage CORAL.  
Lab Door Closed.

## **Problems?**

Send [mtim@mit.edu](mailto:mtim@mit.edu) an email immediately.  
It is OK to try to resolve the problem by restarting the system.

Note: If laser interlock is tripped the system will have to be shutdown and restarted.

## **Shutdown the system:**

Close the shutter in the software.  
Close the software.  
Turn the Main Power Off.  
Turn the Lasers (Key) Off.

## **Restart The System After Power Was Shutoff:**

After a power outage or a tripped laser interlock you will have to restart the system.  
Turn the Main Power ON.  
Open the software (Wire).  
Choose the top option (initialize Motors).  
Turn on the Lasers (Key).

## **Raman Polarization Measurement Configuration Modes**

Polarized –  $\frac{1}{2}$  lambda waveplate - Linear

Depolarized –  $\frac{1}{4}$  lambda waveplate - Circular

### Isotropic Samples

Depolarized Raman

Incident – Depolarized

Analyzer – Polarized X or Y

### Anisotropic Samples

Polarized Raman

Incident – Polarized X and Y

Analyzer – Polarized X and Y

Rotatable Sample Stage

Angle Resolved Polarized Raman

Incident – Linear Polarized ( $\frac{1}{2}$  waveplate) X and Y

Analyzer – Linear Polarized ( $\frac{1}{2}$  Waveplate) 0-190 degrees

Rotatable Sample Stage

**Vendor Contacts:**

<http://www.renishaw.com/en/invia-performance--9452>

**Sample info:**

**Other Instruments:**

Kaiser Hololab 5000 Modular Research Raman Microscope & Spectrometer  
785nm Excitation

Harvard CNS

MIT Institute for Soldier Nanotechnology

**ACKNOWLEDGMENT OF SUPPORT:**

CMSE requires that all publications and patents resulting from work done under this project or with the assistance of the CMSE Facilities, be acknowledged. Please be certain that one (or a combination) of the following statements is included in papers with support from CMSE:

*"This work was supported primarily by the MRSEC Program of the National Science Foundation under award number DMR - 1419807."*

*-or-*

*"This work was supported in part by the MRSEC Program of the National Science Foundation under award number DMR - 1419807."*

*-or-*


"This work made use of the Shared Experimental Facilities supported in part by the MRSEC Program of the National Science Foundation under award number DMR - 1419807."

## **Instrument Specifics:**

### **Renishaw Invia Reflex Raman Confocal Microscope**

This is a fully automated confocal Raman microscope system with three laser excitation lines a motorized stage and open sample compartment to accommodate cryostats. A Temperature controlled microscope stage is available for temperature dependent measurements.

### **Raman Reflex Beampath**



### **Data Acquisition Modes:**

- Point Mapping
- Imaging - Stream HR fast imaging
- Depth Profiling
- Polarization
- Time Dependent
- Temperature Dependent
- Photoluminescence

### **Laser Sources:**

473nm	25mw Linear & Circular
532nm	50mw automated rotate linear
785nm	100mw Linear & Circular

### **Detector:**

1" CCD array, Deep Depletion for enhanced IR Range  
400-1050nm, 1024 X 256 pixel

### **Gratings:**

600l/mm, 1200l/mm, 2400l/mm

### **Rayleigh Filters:**

473nm edge      10%@150cm<sup>-1</sup> shift, for Photoluminescence measurements



532nm eclipse <15cm-1 shift, Stokes & Antistokes  
532nm edge ~100cm-1,  
785nm edge 10%@50cm-1 Rectangle w/o aperture

**Microscope:**

Leica Optical Microscope  
Infinity Corrected, 200mm Tube Length  
Free Space with enclosure

**Objectives:**

Semi par-focal

Objective	Opening	WD	NA
10X	4mm	12mm	0.25
LWD50X	2mm	8mm	0.50
50X	4mm	350um	0.75
LWD100X	2mm	~3mm	0.75
100X	2mm	0.1mm	0.90

Reflected Light Polarization

**Microscope Stages:**

Renishaw HSES motorized stage

Step Size: 100nm

Stage Travel: 112mm x 76mm

**Rotating Microscope stage** - polarization experiments

**Temperature controlled microscope stage**

Linkam THMS600 (-196deg C to 600deg C).

Macro Sampling Set:

Cuvettes, Powders , Capillaries, Etc.

Computer & Monitor

HP Computer

Monitor 22"

Reference Spectra:

Organic Polymers, Inorganic materials, Sample substrates

Features:

Data Analysis Software available to users.

\*\*\*\*\*End\*\*\*\*\*